SolFlex Cast Wax

Gebrauchsinformation

Deutsch

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English

Description:

SolFlex Cast Wax is a light-curing resin for 3D printing of cast patterns with DLP® systems. It can be used for dental technology applications (3D printing of cast or pressable patterns), hearing aid applications (3D printing of cast forms for silicone earmolds) as well as for jewelry (3D printing of cast patterns).

SolFlex Cast Wax is optimized for SolFlex 3D printers (wavelength: 385 nm).

Composition:

Composition of acrylates, initiators and dyes.

Application:

- 3D printing of metal cast patterns in dental technology
- 3D printing of cast form for silicone earmolds
- 3D printing jewelry cast patterns

Processing

1. Manufacturing process

1.1. Resin must be stirred thoroughly with a spatula in the container before pouring into the material reservoir of the machine.

1.2. The manufacturing process has to be performed according to the machine manufacturer's specifications.

1.3. Chosen parameters have to ensure sufficient adhesion between the layers.

2. Cleaning

Observe the references in chapter „Warning notes“ when using isopropanol.

2.1. Conduct cleaning with isopropanol in a separate beaker. Put the beaker in an ultrasonic bath for max. 5 mins. On filigree parts, renouncing the use of an ultrasonic bath is recommended. Clean the objects carefully with isopropanol in a spray bottle and a brush instead.

2.2. Afterwards, blow dry the parts with compressed air and let them dry thoroughly for 30 mins. If available, objects can be stored in a heating cabinet at 37°C.

2.3. If necessary, remove support structure.

2.4. In case of insufficient cleaning (shiny areas on the object after drying) repeat steps 2.1 and 2.2.

3. Post Curing

3.1. Parts have to be dried completely and cooled down to room temperature for post curing.

3.2. When using stroboscope devices (such as Otoflash G171, NK-Optik) cure top and bottom of the object with 2000 flashes each. Devices with inert atmosphere prevent the creation of inhibition layer.

Polymerisation units equipped with two UV-A and two blue light bulbs (such as Dentalux 4, Mega-Physik) require 2 x 10 minutes. Using an alternative light source may affect the final properties of the product. Please notice that the light sources need routine

maintenance following manufacturer's instructions.

4. Final Cleaning

4.1. In order to remove the inhibition layer clean the objects again according to step 2.1.

4.2. Dry objects according to step 2.2.

Warning notes:

Please always pay attention to the current safety data sheet!

Avoid breathing vapours, use extraction unit. Avoid contact respectively wear protective gloves, protective clothing, eye protection. If on skin: wash with plenty of water. If skin irritation occurs: Get medical advice. Uncured material is harmful to environment, avoid release to the environment.

Advice on handling isopropanol during cleaning process:

Isopropanol is a hazardous good. You must use a suitable extraction system when working with isopropanol in order to prevent the creation of an explosive atmosphere. Otherwise, explosion-proof devices must be utilized. Consequently observe the valid safety data sheets as well as country-specific legal regulations when handling hazardous substances.

Storage instructions:

Protect uncured material from light. May polymerize spontaneously. Do not use after the expiration date. Store between 4°C and 25°C.

Disposal instructions:

Dispose of contents/container in accordance with local regulations.

Delivery units:

Art.-Nr. REBFCW02 **SolFlex Cast Wax** 500 g
Art.-Nr. REBFCW01 **SolFlex Cast Wax** 1000 g

The product was developed for use in dentistry and audiology and must be used according to the instructions for use. The manufacturer does not accept any liability for damage caused by improper use of system and material. Furthermore, the user must independently check the material for its suitability and applicability for the intended purposes before using it, especially if he deviates from the process described in the instructions for use, or if his process is not described in the instructions for use.

DeltaMed

Manufactured by:
DeltaMed GmbH
Raiffeisenstraße 8a
D-61169 Friedberg
Phone: +49 6031-7283-0
Fax: +49 6031-7283-29
www.deltamed.de

Distributed by:
W2P Engineering GmbH
Am Campus 1
A-3400 Klosterneuburg
Phone: +43 1 3062 857
www.way2production.at