## SolFlex Cast Red

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## SolFlex Cast Red

Light-curing resin based on (meth)acrylate, for cast parts in dental precision casting, burn-out

for **SolFlex Printer** (DLP, 385 nm)

### Important notes

This is a medical device, only to be used by trained specialist personnel only.

### Processing

- ▶ Phosphate-bonded investment material can be used for fabrication of the casting mould. A pre-heating temperature of approx. 800 °C is to be maintained for approx. 30-45 minutes, depending on the size of the casting mould, regardless of the final casting temperature.
- ▶ The properties of the final product depend, among other things, on post-processing. Therefore it must be ensured that the light unit is in an orderly condition and that the object are completely cured (observe process description on page 2).
- ▶ Processing temperature 23 °C ± 2 °C.

### Safety

- ▶ Please follow the instructions on the safety data sheet!
- ▶ Be sure to use personal protective equipment (protective gloves and protective glasses) during processing.
- ▶ Avoid direct contact with the liquid material and the components prior to post-curing. Irritating to eyes and skin (sensitisation is possible).
- ▶ After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- ▶ After contact with skin wash immediately with water and soap.

### Storage

- ▶ **SolFlex Cast Red** is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.

### General

- ▶ Always keep container tightly sealed, immediately close the container carefully after each use.

### Indication:

Fabrication of residue-free combustible cast objects for precision casting technology employing the DLP process.

**Processing:**  
at 23 °C ± 2 °C

### Storage

### Ordering information:

<b>SolFlex Cast Red</b> 1000 g bottle	<b>03902</b>
<b>SolFlex Model Blue</b> 1.000 g bottle	<b>04143</b>
<b>SolFlex Model</b> 1.000 g bottle	
Sand	<b>03866</b>
Ivory	<b>03867</b>
Grey	<b>03900</b>
<b>SolFlex Prov</b> 1.000 g bottle	
A1	<b>04149</b>
A2	<b>04150</b>
A3	<b>04151</b>
500 g bottle	
A1	<b>04146</b>
A2	<b>04147</b>
A3	<b>04148</b>
<b>SolFlex Tray Green</b> 1.000 g bottle	<b>04145</b>
<b>SolFlex Ortho Clear</b> 1.000 g bottle	<b>04144</b>
<b>SolFlex Med Clear</b> 1.000 g bottle	<b>03901</b>

**SolFlex Cast Red** was developed for use in the dental laboratory and must be used in accordance with the instructions for processing and safety. No liability is accepted for damages caused by faulty or improper use of system and materials.

# SolFlex Cast Red

## Manufacturing process

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

## Construction process

Generation of a Print Job complying with machine and material parameters

## Post-processing

If possible, post-processing should commence immediately following the construction process. After raising the platform, a drip time of approx. 10 minutes is recommended.

## Pre-cleaning

Remove construction components from the platform and clean in a separate vessel with isopropyl alcohol (purity  $\geq 98\%$ ) for 3 min. in an ultrasonic bath.

## Cleaning

Then thoroughly clean the openings, cavities and gap areas, if necessary also with compressed air, and, if applicable, remove the construction components carefully from the support structure.

## Main cleaning process

The main cleaning process is performed in a separate vessel with fresh isopropyl alcohol (purity  $\geq 98\%$ ) for 3 min. in an ultrasonic bath. Prior to drying, check for residues.

## Drying

Heat the construction components for 30 min. in an oven to approx. 40°C to remove the solvent residues from the cleaning process.

## Post-exposure

Post-exposure is performed with a xenon photoflash unit Otofash G171 with 2 x 2000 flashes under inert gas conditions (nitrogen), rotate components in between.