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Produktbeschreibung

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SolFlex Model Beige

Inhaltsstoffe

Tripropylenglycoldiacrylat, Hydroxypropylmethacrylat, Diphenyl(2,4,6-trimethylbenzoyl)phosphinoxid

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
Instructions for use

Product description

SolFlex Model Beige is a light-curing acrylic for the generative production of dental models (also for the thermofforming technique) for DLP printers with a 385 nm UV LED light source.

Indication Dental working and presentation models

Colour Beige

Processing

SolFlex Model Beige has been optimised for the SolFlex 3D printing system. Following the printing process, the printed objects must be cleaned, dried and post-cured with the Otoflash G171 xenon photo-flash unit to guarantee the requisite product characteristics.

Note

Use separate material containers and cleaning baths for each printing material so as to avoid the possibility of cross-contamination.

- Recommended wall thickness if the model is not filled: 3 mm
- For precision work, orient the corresponding models horizontally and cores vertically on the building platform.
- For more efficient printing, deep drawn and antagonist models can be positioned vertically on the building platform.
- SolFlex Model Beige can be printed in the following layer thickness: 50 µm
- We recommend allowing a dripping time of approx. 10 minutes once printing is complete.
- Finish the printed objects without delay so as to avoid distortion of their shape and guarantee the required precision.
- Any excess material from the material container following printing should be returned to its original packaging. Please ensure that excess material is free from contaminants and polymerised components.

POST-PROCESSING

Cleaning

We recommend a solvent-resistant cleaning mat for cleaning. Detach the printed objects carefully from the building platform. Cleaning should be performed with isopropanol (≥ 98 % purity) in an unheated ultrasound bath. The operating instructions from the manufacturers must be observed and followed.

The printed objects should be cleaned in two stages:

- | | |
|---------------------------------|---|
| 1) Pre-cleaning | Pre-clean the printed object for 2 minutes in a reusable isopropanol ultrasound bath. |
| Note | The bath's cleaning efficacy decreases the more it is used. If the cleaning performance deteriorates, the bath should be replaced. |
| 2) Main cleaning process | Then clean the printed objects for a further 2 minutes in a fresh isopropanol ultrasound bath. Should there still be resin residue on the printed object after the main cleaning process, the printed object can be cleaned again for max. 1 minute in the isopropanol ultrasound bath. If no more resin residues are visible on the printed object, the printed object should be dried carefully using compressed air. |

Finishing

Remove and sand the support structures carefully using a rotary instrument and without exerting pressure prior to the post-curing. Use a suction device. Remove any remaining resin dust carefully using compressed air. Then rinse printed objects for a few seconds with fresh isopropanol and remove any dust using a fine, clean brush if necessary. Then dry the printed objects carefully using compressed air.

Post-exposure

The post-exposure is performed with an Otoflash G171 xenon photoflash unit with 2 x 2,000 flashes. No protective gas atmosphere is required. However, ensure that the printed objects do not cast shadows over or touch each other, as this could result in incomplete post-curing. Following the first 2,000 flashes, allow a cooling-off period of at least 2 minutes with the lid open. Insufficient cooling can distort the shape of the printed object. Once the cooling-off period has elapsed, the printed objects should be turned and light-cured with a further 2,000 flashes.

Cleaning of printed objects

Any excess foreign materials (e.g., dental wax) can be removed by steaming them off the printed objects. Avoid excessive heating. For example, do not place printed objects in the dental wax extraction unit, as this can distort their shape.

Safety notice

- Please follow the instructions of SolFlex Model Beige safety data sheet!
- SolFlex Model Beige contains (meth)acrylates. SolFlex Model Beige should not be used in case of known hypersensitivities (allergies) to any of these ingredients. SolFlex Model Beige may only be used in its completely cured state. Pay attention to the finishing process.

Storage

Store at temperatures of 15 °C to 28 °C. Seal bottle again immediately after use. The material will cure if exposed to light. Do not use after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and / or written information has been given to the best of our knowledge and without obligation. Our information and / or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

SolFlex Model Beige

Contains

Tripropyleneglycol Diacrylate, Hydroxypropyl Methacrylate, diphenyl (2,4,6-trimethylbenzoyl)phosphine oxide


Hazard statements

Causes skin irritation. Causes serious eye irritation. May cause an allergic skin reaction. Harmful to aquatic life with long lasting effects.

Precautionary statements

Avoid breathing spray. Avoid release to the environment. Wear protective gloves / eye protection / face protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Specific treatment (see in this instructions for use). Dispose of contents / container in accordance with local / regional / national / international regulations.

1000 g

Manufactured by

VOCO GmbH
Anton-Flettner-Str. 1-3
27472 Cuxhaven · Germany
Phone +49 (4721) 719-0
www.voco.dental

Distributed by

W2P Engineering GmbH
Am Campus 1
3400 Klosterneuburg · Austria
Phone +43 (1) 306 28 57
www.way2production.at

VC 60 007121 E1 0618 99