SolFlex Model Sand, Grey, Ivory

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SolFlex Model Sand, Grey, Ivory

Light-curing resin based on (meth)acrylate, for the generative fabrication of dental models

for **SolFlex** Printer (DLP, 385 nm)

Important notes

This is a medical device, only to be used by trained specialist personnel only.

Processing

- ▶ The properties of the final product depend, among other things, on post-processing. Therefore it must be ensured that the light unit is in an orderly condition and that the models are completely cured (observe process description on page 2).
- ▶ Processing temperature 23 °C ± 2 °C.

Safety

- ▶ Please follow the instructions on the safety data sheet!
- ▶ Be sure to use personal protective equipment (protective gloves and protective glasses) during processing.
- ▶ Avoid direct contact with the liquid material and the components prior to post-curing. Irritating to eyes and skin (sensitisation is possible).
- ▶ After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- ▶ After contact with skin wash immediately with water and soap.

Storage

- ▶ **SolFlex Model** is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.

General

- ▶ Always keep container tightly sealed, immediately close the container carefully after each use.


Indication:

3D Print working and presentation models for the entire spectrum of models in dentistry

Processing:

at 23 °C ± 2 °C

Storage

Ordering information:

SolFlex Model	
1.000 g bottle	
Sand	03866
Ivory	03867
Grey	03900
SolFlex Model Blue	
1.000 g bottle	04143
SolFlex Prov	
1.000 g bottle	
A1	04149
A2	04150
A3	04151
500 g bottle	
A1	04146
A2	04147
A3	04148
SolFlex Tray Green	
1.000 g bottle	04145
SolFlex Ortho Clear	
1.000 g bottle	04144
SolFlex Med Clear	
1.000 g bottle	03901
SolFlex Cast Red	
1000 g bottle	03902

SolFlex Model was developed for use in the dental laboratory and must be used in accordance with the instructions for processing and safety. No liability is accepted for damages caused by faulty or improper use of system and materials.

SolFlex Model Sand, Grey, Ivory

Manufacturing process

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

Construction process

Generation of a Print Job complying with machine and material parameters

Post-processing

If possible, post-processing should commence immediately following the construction process. After raising the platform, a drip time of approx. 10 minutes is recommended.

Pre-cleaning

Remove construction components from the platform and clean in a separate vessel with isopropyl alcohol (purity $\geq 98\%$) for 3 min. in an ultrasonic bath.

Cleaning

Then thoroughly clean the openings, cavities and gap areas, if necessary also with compressed air, and, if applicable, remove the construction components carefully from the support structure.

Main cleaning process

The main cleaning process is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq 98\%$) for 3 min. in an ultrasonic bath. Prior to drying, check the openings and additional bore holes for residues.

Drying

Heat the construction components for 30 min. in an oven to approx. 40°C to remove the solvent residues from the cleaning process.

Post-exposure

Post-exposure is performed with a xenon photoflash unit Otoflash G171 with 2 x 2000 flashes under inert gas conditions (nitrogen), rotate components in between.