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Produktbeschreibung

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SolFlex Surgical Guide

Inhaltsstoffe

Diphenyl(2,4,6,-trimethylbenzoyl)phosphinoxid

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Instructions for use

Product description

SolFlex Surgical Guide is a light-curing acrylic for the generative production of biocompatible dental drilling templates (medical device class IIa) for DLP printers with a 385 nm UV LED light source. Drilling templates made of SolFlex Surgical Guide can be steam sterilised at 134 °C for max. 5 minutes.

Indication Dental drilling templates

Colour: clear

Processing

SolFlex Surgical Guide has been optimised for the SolFlex 3D printing system. Following the printing process, the printed objects must be cleaned, dried and post-exposed with the Otofash G171 xenon photoflash unit to guarantee the requisite product characteristics.

Note

Use separate material containers and cleaning baths for each printing material so as to rule out the possibility of cross-contamination.

- Minimum wall thickness for construction: 3 mm
- SolFlex Surgical Guide can be printed in the following layer thickness: 50 µm
- We recommend allowing a dripping time of approx. 10 minutes once printing is complete.
- Finish the printed objects without delay so as to avoid distortion of their shape and guarantee the required precision.
- Any excess material from the material container following printing should be returned to its original packaging. Please ensure that excess material is free from contaminants and polymerised components.

POST-PROCESSING

Cleaning

We recommend a solvent-resistant cleaning mat for cleaning. Detach the printed objects carefully from the build platform. Cleaning should be performed with isopropanol (≥ 98 % purity) in an unheated ultrasound bath. The operating instructions from the manufacturers must be observed and followed.

The printed objects should be cleaned in two stages.

- 1) Pre-cleaning** Pre-clean the printed object for 2 minutes in a reusable isopropanol ultrasound bath.
Note The bath's cleaning efficacy decreases the more it is used. If the cleaning performance deteriorates, the bath should be replaced.
- 2) Main cleaning process** Then clean the printed objects for a further 2 minutes in a fresh isopropanol ultrasound bath. Should there still be resin residue on the printed object after the main cleaning process, the printed object can be cleaned again for max. 1 minute in the isopropanol ultrasound bath. If no more resin residues are visible on the printed object, the printed object should be dried carefully using compressed air.

Finishing

Remove and sand the support structures carefully using a rotary instrument and without exerting pressure prior to the post-exposure. Use a suction device. Remove any remaining resin dust carefully using compressed air. Then rinse printed objects for a few seconds with fresh isopropanol and remove any dust using a fine, clean brush if necessary. Then dry the printed objects again carefully using compressed air.

Bonding of drill sleeves

Coat the exterior of the sleeve and the intended hole in the drilling template with a thin layer of liquid SolFlex Surgical Guide prior to the post-exposure. Insert the sleeves and remove any excess material to ensure marginal integrity. Then post-expose the drilling template as described under "Post-exposure".

Post-exposure

The post-exposure is performed with an Otofash G171 xenon photoflash unit with 2 x 2,000 flashes. No protective gas atmosphere is required. However, ensure that the printed objects do not cast shadows over or touch each other, as this could result in incomplete post-curing.

Following the first 2,000 flashes, allow a cooling-off period of at least 2 minutes with the lid open. Insufficient cooling can distort the shape of the printed object. Once the cooling-off period has elapsed, the printed objects should be turned and light-cured with a further 2,000 flashes.

Sterilisation

The sterilisation process is validated according to DIN EN ISO 17664 and can be used with the equipment that complies within EN 13060. The instructions for use from the manufacturer have to be observed and followed. The user is responsible for the sterility of the product. Only use completely post-exposed and cooled drilling templates to undergo the steam sterilisation! Otherwise a deformation of the printed object cannot be excluded.

Technical data of the sterilisation process

134 °C	2,07 - 2,17 bar	max. 5 min
--------	-----------------	------------

Disinfection

After use, drilling templates produced from SolFlex Surgical Guide can be disinfected with alcohol- or aldehyde-based disinfectant solutions (e.g., ethanol, MD 520 from DÜrr, Cavex Impre Safe from Cavex). Observe the manufacturer's instructions for use.

Information, precautionary measures

SolFlex Surgical Guide may only be used intraorally in its completely cured state. Pay attention to the finishing process. SolFlex Surgical Guide contains methacrylates. SolFlex Surgical Guide should not be used in case of known hypersensitivities (allergies) to any of these ingredients.

Storage

Store between 15 °C and 28 °C. Seal bottle again immediately after use. The material will cure if exposed to light. Do not use after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

SolFlex Surgical Guide

Contains

Diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide

Hazard statements

May cause an allergic skin reaction. Harmful to aquatic life with long lasting effects.

Precautionary statements

Avoid breathing spray. Avoid release to the environment. Wear protective gloves. Specific treatment (see in this instructions for use). If skin irritation or rash occurs: Get medical advice / attention. Dispose of contents / container in accordance with local / regional / national / international regulations.

1000 g

Achtung
Warning

Manufactured by

VOCO GmbH
Anton-Flettner-Str. 1-3
27472 Cuxhaven · Germany
Phone +49 (4721) 719-0
www.voco.dental

Distributed by

W2P Engineering GmbH
AM Campus 1
3400 Klosterneuburg · Austria
Phone +43 (1) 306 28 57
www.way2production.at

VC 60 007122 E1 0618 99