SolFlex Tray Green

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SolFlex Tray Green

Light-curing resin based on (meth)acrylate, biocompatible, for the generative fabrication of individual and functional impression trays for **SolFlex Printer** (DLP, 385 nm)

Indication:

Individual impression trays, functional impression trays and base plates

Important notes

This is a medical device, only to be used by trained specialist personnel.

Processing

- ▶ The properties of the final product depend, among other things, on post-processing. Correct post-exposure is important for biocompatibility. Therefore it must be ensured that the light unit is in an orderly condition and that the moulds are completely cured (observe process description on page 2).
- ▶ Minimal material thickness of the tray 2.5 mm, in the grip area 5.0 mm
Minimum material thickness of base plates 2.5 mm
- ▶ Maximum curing depth* at direct post-exposure:
clear: 6 mm
*In case of large objects and exposure on both sides, the material thickness can be up to 12 mm (Example **SolFlex Tray Green** – with a curing depth of 6 mm).
Do not use heat-based methods for disinfection or sterilisation. This could possibly deform the work-piece.
- ▶ Polish surface mechanically.
- ▶ Processing temperature 23 °C ± 2 °C

Safety

- ▶ Please follow the instructions on the safety data sheet!
- ▶ Be sure to use personal protective equipment (protective gloves and protective glasses) during processing.
- ▶ Avoid direct contact with the liquid material and the components prior to post-curing. Irritating to eyes and skin (sensitisation is possible).
- ▶ After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- ▶ After contact with skin wash immediately with water and soap.
- ▶ Biocompatibility is only guaranteed with complete polymerisation.

Storage

- ▶ **SolFlex Tray Green** is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.

General

- ▶ Always keep container tightly sealed, immediately close the container carefully after each use.

Contraindication

Contains (meth)acrylics and phosphine oxide.

Some ingredients of **SolFlex Tray Green** may cause allergic reactions in predisposed persons. In such cases refrain from using the product. **SolFlex Tray Green** only insert intraorally in completely polymerised state.

Adverse effects

Product may cause allergic reactions.

Processing:

at 23 °C ± 2 °C

Storage

Ordering information:

| | |
|--|--------------|
| SolFlex Tray Green 1.000 g bottle | 04145 |
| SolFlex Model Blue 1.000 g bottle | 04143 |
| SolFlex Model 1.000 g bottle | |
| Sand | 03866 |
| Ivory | 03867 |
| Grey | 03900 |
| SolFlex Prov 1.000 g bottle | |
| A1 | 04149 |
| A2 | 04150 |
| A3 | 04151 |
| 500 g bottle | |
| A1 | 04146 |
| A2 | 04147 |
| A3 | 04148 |
| SolFlex Ortho Clear 1.000 g bottle | 04144 |
| SolFlex Med Clear 1.000 g bottle | 03901 |
| SolFlex Cast Red 1000 g bottle | 03902 |

SolFlex Tray Green was developed for use in dental laboratory and must be used in accordance with the instructions for processing and safety. DETAX will not be responsible for damages caused by faulty or improper use of system and materials.

Made in Germany

08/2017

SolFlex Tray Green

Manufacturing process

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

Construction process

Generation of a Print Job complying with machine and material parameters

Post-processing

If possible, post-processing should commence immediately following the construction process. After raising the platform, a drip time of approx. 10 minutes is recommended.

Pre-cleaning

Remove construction components from the platform and clean in a separate vessel with isopropyl alcohol (purity $\geq 98\%$) for 3 min. in an ultrasonic bath.

Cleaning

Then thoroughly clean the printing objects, if necessary also with compressed air, and, if applicable, remove the construction components carefully from the support structure.

Main cleaning process

The main cleaning process is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq 98\%$) for 3 min. in an ultrasonic bath. Prior to post-exposure, check the openings and additional bore holes for residues. Then blow off with compressed air.

Post-exposure

Post-exposure is performed with a xenon photoflash unit (e.g. Otofash G171) with 2 x 2000 flashes under inert gas conditions (nitrogen), rotate components in between.

Surface processing

The surface can be polished as required.